

Lateral Flow Assay Development Kit (Starter Kit)



LFAD10007			LFAD10008		
Items	qty	remark	Items	qty	remark
NC membrane	1	100m/roll, 2.5cm, 120	NC membrae	2	100m/roll, 2.5cm, 90-120,
sample pad	100	16*300mm	sample pad	100	16*300mm
conjugate pad	100	10*300mm	conjugate pad	100	10*300mm
absorbent pad	100	18*300mm	absorbent pad	100	18*300mm
backing card	100	300*60mm (20+25+15)	backing card	100	300*60mm (20+25+15)
blood separator	100	10*30mm	blood separator	100	10*30mm
lateral flow cassette	100	20*70mm; 3mm	lateral flow cassette	100	20*70mm; 3mm
foil bag	200	125*65mm; blank	foil bag	200	125*65mm; blank

LFAD10009			LFAD10010		
Items	qty	remark	Items	qty	remark
NC membrane	3	100m/roll, 2.5cm, 90-120-140	NC membrae	4	100m/roll, 2.5cm, 90-120-140-180
sample pad	100	16*300mm	sample pad	100	16*300mm
conjugate pad	100	10*300mm	conjugate pad	100	10*300mm
absorbent pad	100	18*300mm	absorbent pad	100	18*300mm
backing card	100	300*60mm (20+25+15)	backing card	100	300*60mm (20+25+15)
blood separator	100	10*30mm	blood separator	100	10*30mm
lateral flow cassette	100	20*70mm; 3mm	lateral flow cassette	100	20*70mm; 3mm
foil bag	200	125*65mm; blank	foil bag	200	125*65mm; blank