



Lateral Flow Backing Card

For lateral flow test strip manufacturing

Daily produce 500,000 sheets / worldwide shipment



Structure



Split line illustration



Parameter

Item	Lateral Flow Backing Card				
Peeling (thickness)	0.1mm				
Adhesive (thickness)	0.03mm	0.03mm	0.03mm	0.03mm	0.03mm
PVC (thickness)	0.27mm	0.32mm	0.37mm	0.42mm	0.47mm
Total thickness	0.3mm	0.35mm	0.4mm	0.45mm	0.5mm
Width	Customized				
Split line					

Remark: total thickness exclude peeling thickness

Commonly used dimension

Sheet format	Roll format
300*60mm	60mm*100M
300*70mm	70mm*100M
300*80mm	80mm*100M
Customized	

Commonly used dimension

No.	Dimension	Format	Thickness (mm)	Application
CPVC0001	300*60mm	Sheet	0.3	Colloidal gold immunochromatography
CPVC0002	300*60mm		0.35	
CPVC0003	300*60mm		0.4	
CPVC0004	300*60mm		0.45	
CPVC0005	300*60mm		0.5	
CPVC0006	300*70mm		0.3	
CPVC0007	300*70mm		0.35	
CPVC0008	300*70mm		0.4	
CPVC0009	300*70mm		0.45	
CPVC0010	300*70mm		0.5	
CPVC0011	300*80mm		0.3	
CPVC0012	300*80mm		0.35	
CPVC0013	300*80mm		0.4	
CPVC0014	300*80mm		0.45	
CPVC0015	300*80mm		0.5	

No.	Dimension	Format	Thickness (mm)	Application
CPVC0016	60mm*100M	Roll	0.3	Colloidal gold immunochromatography
CPVC0017	60mm*100M		0.35	
CPVC0018	60mm*100M		0.4	
CPVC0019	60mm*100M		0.45	
CPVC0020	60mm*100M		0.5	
CPVC0021	70mm*100M		0.3	
CPVC0022	70mm*100M		0.35	
CPVC0023	70mm*100M		0.4	
CPVC0024	70mm*100M		0.45	
CPVC0025	80mm*100M		0.5	
CPVC0026	80mm*100M		0.3	
CPVC0027	80mm*100M		0.35	
CPVC0028	80mm*100M		0.4	
CPVC0029	80mm*100M		0.45	
CPVC0030	80mm*100M		0.5	

No.	Dimension	Format	Thickness (mm)	Application
FPVC0001	300*60mm	Sheet	0.3	Fluorescence immunochromatography
FPVC0002	300*60mm		0.35	
FPVC0003	300*60mm		0.4	
FPVC0004	300*60mm		0.45	
FPVC0005	300*60mm		0.5	
FPVC0006	300*70mm		0.3	
FPVC0007	300*70mm		0.35	
FPVC0008	300*70mm		0.4	
FPVC0009	300*70mm		0.45	
FPVC0010	300*70mm		0.5	
FPVC0011	300*80mm		0.3	
FPVC0012	300*80mm		0.35	
FPVC0013	300*80mm		0.4	
FPVC0014	300*80mm		0.45	
FPVC0015	300*80mm		0.5	

No.	Dimension	Format	Thickness (mm)	Application
FPVC0016	60mm*100M	Roll	0.3	Fluorescence immunochromatography
FPVC0017	60mm*100M		0.35	
FPVC0018	60mm*100M		0.4	
FPVC0019	60mm*100M		0.45	
FPVC0020	60mm*100M		0.5	
FPVC0021	70mm*100M		0.3	
FPVC0022	70mm*100M		0.35	
FPVC0023	70mm*100M		0.4	
FPVC0024	70mm*100M		0.45	
FPVC0025	80mm*100M		0.5	
FPVC0026	80mm*100M		0.3	
FPVC0027	80mm*100M		0.35	
FPVC0028	80mm*100M		0.4	
FPVC0029	80mm*100M		0.45	
FPVC0030	80mm*100M		0.5	