

IVD Consumables for Lateral Flow Assay Development and Manufacturing

Product	Material	Specification	Order No.	Application	Remark
Lateral Flow Assay Consumables	Glass Fiber Membranes	300x200mm	S-100	Sample Pad	Supply by sheet or by roll
		300x200mm	S-80	Sample Pad	
		300x200mm	SC-80	Sample & Conjugate Pad	
		300x200mm	SC-75	Sample & Conjugate Pad	
		300x200mm	C-55	Conjugate Pad	
		300x200mm	C-40	Conjugate Pad	
	Blood Separation Membranes	300x71mm	BLM-55	blood separation	
	Polyester Fiber Membranes	300x200mm	PM-80	Conjugate Pad	
		300x200mm	PM-90	Conjugate Pad	
	Wood pulp/Cotton pulp	300x200mm	AP-F270	Absorbent Pad	
		300x200mm	AP-S270	Absorbent Pad	
		300x200mm	AP-D370	Absorbent Pad	
Lateral Flow Backing Card	PVC	60 x 300mm	PVC036003	Thickness 0.3mm	Supply by sheet or by roll; For lateral flow assay test development and manufacturing. Use as backing on rapid test, suitable for Qualitative and Quantitative IVD diagnostics test kit.
		60 x 300mm	PVC036004	Thickness 0.4mm	
		60 x 300mm	PVC036005	Thickness 0.5mm	
		60 x 300mm	PVC036006	Thickness 0.6mm	
		70 x 300mm	PVC037003	Thickness 0.3mm	
		70 x 300mm	PVC037004	Thickness 0.4mm	
		70 x 300mm	PVC037005	Thickness 0.5mm	
		70 x 300mm	PVC037006	Thickness 0.6mm	
		80 x 300mm	PVC038003	Thickness 0.3mm	
		80 x 300mm	PVC038004	Thickness 0.4mm	
		80 x 300mm	PVC038005	Thickness 0.5mm	
		80 x 300mm	PVC038006	Thickness 0.6mm	
		60 x 300mm	PVC036045	Thickness 0.45mm	
		70 x 300mm	PVC037045	Thickness 0.45mm	
		80 x 300mm	PVC038045	Thickness 0.45mm	
		60 x 300mm	PVC036056	Thickness 0.56mm	
		70 x 300mm	PVC037056	Thickness 0.56mm	
		80 x 300mm	PVC038056	Thickness 0.56mm	
	PS	60 x 300mm	PS0360029	Thickness 0.29mm	
		60 x 300mm	PS0360041	Thickness 0.41mm	
		70 x 300mm	PS0370029	Thickness 0.29mm	
		70 x 300mm	PS0370041	Thickness 0.41mm	
		80 x 300mm	PS0380029	Thickness 0.29mm	
		80 x 300mm	PS0380041	Thickness 0.41mm	

IVD Consumables for Lateral Flow Assay Development and Manufacturing

	PET	60 x 300mm	PET036003	Thickness 0.3mm	
		60 x 300mm	PET036035	Thickness 0.35mm	
		70 x 300mm	PET037003	Thickness 0.3mm	
		70 x 300mm	PET037035	Thickness 0.35mm	
		80 x 300mm	PET038003	Thickness 0.3mm	
		80 x 300mm	PET038035	Thickness 0.35mm	
Nitrocellulose membranes	NC	2.0cm x 100m	NC90-20	<u>Immunochromatography</u> , <u>Lateral Flow Assay</u>	By roll format
		2.5cm x 100m	NC90-25		
		2.0cm x 100m	NC120-20		
		2.5cm x 100m	NC120-25		
		2.0cm x 100m	NC140-20		
		2.5cm x 100m	NC140-25		
		2.0cm x 100m	NC180-20		
		2.5cm x 100m	NC180-25		
Cassette Case	Plastic	20 x 70mm	CC001-360	strip: 3mm x 60mm	Cover & Bottom
		20 x 70mm	CC002-460	strip: 4mm x 60mm	
Nasopharyngeal Swab and Throat Swab	nylon	3 x 150mm	BSWNP80	Nasopharyngeal Swab	For influenza test, covid testing, hand-foot and mouth test, etc.
		4 x 150mm	BSWOP30	Throat Swab	
viral transport medium	dnase/rnase	10ml	VTM1030M	3mL liquid (viral transport media)	For storage virus of SARS-CoV-2 to covid testing pcr detection or SARS-CoV-2 antigen and antibody rapid test
Foil Pouch	PET12/AL7/PP55	120 x 65mm	FPB12065	thickness 7.5cc	For lateral flow rapid test packaging, SARS-CoV-2 rapid test, Influenza Rapid Test, HCG rapid test, Food Safety Rapid test, etc.
		130 x 80mm	FPB13080		
		140 x 65mm	FPB14065		
		140 x 80mm	FPB14080		