







Test Strip Assembly Machine

- 10 Strips or 5 Strips/Sheet OR Single Strip
- 2024 Innovative Rapid Test Packaging Method
- Save up to **60%** packaging cost



Reduce plastic usage, protect environment

Production Line




4,800mm





4,500mm

4,500mm

1,800mm



Advantages

Innovative Packaging	No Need Items Reduce plastic usage, protect environment	
		No need desiccant
		No need plastic cassette
		No need aluminum foil bag

Parameters

Model	Auto / Semi Auto	Workers	Remark	Capacity
TKAM-502	Fully auto	1	In production line	3000 strips/hour
TKAM-506	Semi auto	6	6 individual machines, strip cutter not integrated	1000 strips/hour


Some workstations




Backing film feeder




Cover film cutter



Strip Cutter



CCD Inspection



Robots transfer strips



Strips stick on backing film

Video: <https://youtu.be/ELH2TCohBe0>

Innovative Test Strip Packaging System

Contact info@antiteck.com

Add: A1-519, XingGang GuoJi, Yingbin Road, Huadu, Guangzhou, China

WhatsApp/WeChat/Mob: +86 133 1286 3972 Web: <https://antiteck.com>